

# Disposable Trocar

## 产品特点 FEATURE

**Universal Sealing** design increases resistance of instrument insertion without using the reducer.

**Threaded Design.** accepts both conventional puncture and spiral puncture. The controlled insertion speed and depth of trocar ensures stable fixation and avoids damage to abdomen.

**Pyramid Tip** protects inter organs when puncturing.



**Large Oriented Aiming Conical Surface** facilitates one-handed operation and instrument insertion into the working channel and easy access to large tissues

**Twist Design**

**Translucent Design** for easier observation of instrument insertion/removal as well as extraction of specimen.

**Additional Positioner** to ensure trocar placement .

## Disposable Trocar

Tip	Cannula	Obturator	Seal Style	Length (mm)	Diameter (mm)
N Pyramid	P Threaded T Screw	V Optical C Non-Optical	S Single Sheet M 4 Sheets	100 150(for Bariatric Port )	1-5 1-10 1-12

No.	REF NO.	No.	REF NO.
1	NPVS -100 -1-5	6	NPCS-100-1-5
2	NPVS -100 -1-10	7	NPCS-100-1-10
3	NPVS-100-1-12	8	NPCS-100-1-12
4	NPVM-100-1-10	9	NPCM-100-1-10
5	NPVM-100-1-12	10	NPCM-100-1-12

No.	REF NO.	No.	REF NO.
1	NPCS -150 -1-5	4	NPVS-150-1-5
2	NPCS -150 -1-10	5	NPVS-150-1-10
3	NPCS-100-1-12		