



DISPOSABLE ANORECTAL STAPLER




Website QR Code

Surgaid 施爱德®
Surgaid, an excellent aid for the surgeons

Heal Haemorrhoids, feel No More Fears!


Product features



Type I



Type II



Specifications

NO	Ref No.	Description
1	WH BP-I-32	PPH, Diameter-32 mm
2	WH BP-I-34	PPH, Diameter-34 mm
3	WH BP-I-36	PPH, Diameter-36 mm
4	WH BP-II-32	TST, Diameter-32 mm
5	WH BP-II-34	TST, Diameter-34 mm
6	WH BP-II-36	TST, Diameter-36 mm

Note: For more product information, please consult our marketing department.



TopGun®

DISPOSABLE CIRCUMCISION STAPLER




Website QR Code

Surgaid 施爱德®
Surgaid, an excellent aid for the surgeons

Use my gun to repair your gun !

Product features



Specifications

NO	Ref No.	Description
1	IIB 26	Diameter-26 mm
2	IIB 29	Diameter-29 mm
3	IIB 34	Diameter-34 mm

Note: For more product information, please consult our marketing department.



DISPOSABLE INTRALUMINAL STAPLER




Website QR Code

Surgaid 施爱德®
Surgaid, an excellent aid for the surgeons

Disposable intraluminal stapler, the safety you can feel!

Product features




Specifications

NO	Ref No.	Description
1	WRB-a 21	Diameter-21 mm
2	WRB-a 26	Diameter-26 mm
3	WRB-a 29	Diameter-29 mm
4	WRB-a 34	Diameter-34 mm

Note: For more product information, please consult our marketing department.



DISPOSABLE LINEAR CUTTER STAPLER AND RELOADS




Website QR Code

Surgaid 施爱德®
Surgaid, an excellent aid for the surgeons

Single Use Sterile Specimen Pouch, the safety you can feel!


Product features



Specifications

Product Code	Staple Length (mm)	Stapler Height (mm)
WKD	60、80	B-3.8(BLUE)、G-4.7(GREEN)

No.	Ref no.	Description
1	WKD-60B	Length-60mm, stapler height-3.8mm
2	WKD-60G	Length-60mm, stapler height-4.7mm
3	WKD-80B	Length-80mm, stapler height-3.8mm
4	WKD-80G	Length-80mm, stapler height-4.7mm

Reload	Staple Length (mm)	Stapler Height (mm)
	60、80	B-3.8(BLUE)、G-4.7(GREEN)

No.	Ref no.	Description
1	ZKD-60B	Length-60mm, stapler height-3.8mm
2	ZKD-60G	Length-60mm, stapler height-4.7mm
3	ZKD-80B	Length-80mm, stapler height-3.8mm
4	ZKD-80G	Length-80mm, stapler height-4.7mm

Note: For more product information, please consult our marketing department.



DISPOSABLE LAPAROSCOPIC LINEAR STAPLER I



Website QR Code

Surgaid 施爱德®

Surgaid, an excellent aid for the surgeons

Disposable laparoscopic linear stapler, the safety you can feel!

Product features

The special "droplets" shaping groove design can ensure B-shape stapler staple formation and provide better hemostatic effect and air tightness effect for clinical operation.

Indicate the direction of the blade walking to ensure safe operation without any accident.

Reset the button manually for firing.

3+1 firing model, manual firing back after 3 step of firing, safer and more stable.

Adaptive joint head can swing within 90°. The joint moves flexible and can be operated by single hand conveniently.

Specifications

Product Code	Device Length	Staple Length	Staple Height
WQD	S - Short (270mm) E - Standard (330mm) L - Long (430mm)	45、60、75	A - 2.5 Y - 3.8 C - 3.5 G - 4.1

No.	REF NO.	Device Length	No.	REF NO.	Device Length	No.	REF NO.	Device Length
1	WQDbS-60A	Short	5	WQDbE-60A	Standard	9	WQDbL-60A	Long
2	WQDbS-60C		6	WQDbE-60C		10	WQDbL-60C	
3	WQDbS-60Y		7	WQDbE-60Y		11	WQDbL-60Y	
4	WQDbS-60G		8	WQDbE-60G		12	WQDbL-60G	

Reload I	Staple Length (mm)	Staple Heights
	45、60、75	A - 2.5 (White) C - 3.5 (Blue) Y - 3.8 (Gold) G - 4.1 (Green)

No.	REF NO.	Staple Length	No.	REF NO.	Staple Length	No.	REF NO.	Staple Length
1	ZQ Nb-45A	45	5	ZQ Nb-60A	60	9	ZQ Nb-75A	75
2	ZQ Nb-45C		6	ZQ Nb-60C		10	ZQ Nb-75C	
3	ZQ Nb-45Y		7	ZQ Nb-60Y		11	ZQ Nb-75Y	
4	ZQ Nb-45G		8	ZQ Nb-60G		12	ZQ Nb-75G	

Note: For more product information, please consult our marketing department.



DISPOSABLE LAPAROSCOPIC LINEAR STAPLER II




Website QR Code

Surgaid 施爱德®
Surgaid, an excellent aid for the surgeons

Disposable laparoscopic linear stapler, the safety you can feel!

Product features



Specifications

Device Code		Device Length	
WQ		S - Short (60mm) E - Standard (160mm) L - Long (260mm)	
No.	REF NO.		
1、2、3	WQF-S、WQF-E、WQF-L		
Reload II	Staple Length	Open Staple Heights	
	30、45 60、70	A - 2.5 (White) C - 3.5 (Blue) G - 4.1 (Gold) V - 4.8 (Green)	
No.	REF NO.	No.	REF NO.
1	ZQ I Ud-45A	5	ZQ I Ud-70A
2	ZQ I Ud-45C	6	ZQ I Ud-70C
3	ZQ I Ud-45G	7	ZQ I Ud-70G
4	ZQ I Ud-45V	8	ZQ I Ud-70V

Note: For more product information, please consult our marketing department.