



# DISPOSABLE TROCAR KIT



Website QR Code

**Surgaid** 施爱德®


*Surgaid, an excellent aid for the surgeons*

**Disposable laparoscopic trocar kit, the safety you can feel!**

## Product features




**Appropriate Sealing**  
offers good air tightness.



**Threaded Design** control  
insertion speed and depth.



**Visualization tip** allows  
excellent identification.



**Smooth conical aiming**  
facilitates easier insertion.



**Simple seal detachment**  
for versatility during specimen  
removal.



**Fixed button** with advanced  
fixation keep scope from moving  
during procedures.

## Specifications

No.	Ref No.	Description
1	NPVS-100-1-5	Optical, 5*100mm
2	NPVS-100-1-10	Optical, 10*100mm
3	NPVS-100-1-12	Optical, 12*100mm
4	NPVS-100-2-B	Optical, Length-100mm (1 Trocar 10mm + 1 Cannula 10mm)
5	NPVS-100-3-A	Optical, Length-100mm (1 Trocar 5mm + 1Trocar 10mm + 1 Cannula 10mm)
6	NPVS-100-4-B	Optical, Length-100mm (1 Trocar 5mm + 1 Trocar 10mm + 1 Trocar 12mm + 1 Cannula 5mm)
7	NPVS-100-4-C	Optical, Length-100mm (1 Trocar 5mm + 1 Cannula 5mm + 1 Trocar 12mm + 1 cannula 12mm)

Note: For more product information, please consult our marketing department.



# DISPOSABLE TROCAR KIT



**Disposable laparoscopic trocar kit, the safety you can feel!**

## Product features




**Universal Sealing**  
increases trocar insertion resistance without using reducer.




**Threaded Design** enables to control puncture speed and depth and fixed firmly that avoids damage to abdomen.



**Pyramid Tip** protects inter organs when puncturing.



**Large aiming tapered surface** facilitates one-handed operation and easy insert working channel.



**Translucent Design** for easier observation of instrument insertion/removal and extraction of specimen.

## Specifications

No.	Ref No.	Description
1	NPCS-100-1-5	Bladeless, 5*100mm
2	NPCS-100-1-10	Bladeless, 10*100mm
3	NPCS-100-1-12	Bladeless, 12*100mm
4	NPCS-100-1-15	Bladeless, 15*100mm
5	NPCS-150-1-5	Bladeless, 5*150mm, For Bariatric Surgery
6	NPCS-150-1-10	Bladeless, 10*150mm, For Bariatric Surgery
7	NPCS-150-1-12	Bladeless, 12*150mm, For Bariatric Surgery

Note: For more product information, please consult our marketing department.



Website QR Code

**Surgaid** 施爱德®

*Surgaid, an excellent aid for the surgeons*



## DISPOSABLE PEDIATRIC TROCAR




Website QR Code

**Surgaid** 施爱德®  
Surgaid, an excellent aid for the surgeons

**Disposable pediatric trocar, the safety you can feel!**

### Minimal invasive




### Specifications

No.	Ref No.	Description
1	IIC-50-1-3	Bladeless, 3*50mm, plastic obturator
2	IIC-70-1-3	Bladeless, 3*70mm, metal obturator
3	IIC-70-1-5	Bladeless, 5*70mm, plastic obturator
4	IIC-100-1-5	Bladeless, 5*100mm, plastic obturator
5	IIV-100-1-5	Bladeless, 5*100mm, metal obturator


Note: For more product information, please consult our marketing department.

## DISPOSABLE THYROIDECTOMY TROCAR



### ■ Indications

Benign diseases such as simple goiter, nodular goiter (or with cystic hyperplasia), or thyroid adenoma.



### ■ Specifications

No.	Ref no.	Description
1	IC-130-1-5	Length-130mm, Diameter-5mm
2	IC-130-1-12	Length-130mm, Diameter-12mm
3	IC-130-3-F	Length-130mm, Diameter(5mm/5mm/12mm)

Note: For more product information, please consult our marketing department.




Website QR Code

**Surgaid** 施爱德®  
Surgaid, an excellent aid for the surgeons



## DISPOSABLE SMOKE FILTER




Website QR Code

**Surgaid** 施爱德®  
Surgaid, an excellent aid for the surgeons

*Breathe healthy and enjoy safety all the time*

### Effective Filter

Contains activated carbon to absorb harmful gases.




### Universal Connector

High compatibility to connect to the trocar.




### Moisture Management

Built-in fluid trap to manage moisture throughout entire procedure.



### Adjustable flow rate

0: Closing velocity    1: 50% velocity  
2: 90% velocity    3: 100% velocity



### Specifications

No.	Ref No.	Description
1	T	Smoke Filter
2	IIVT-100-3-A	Optical, Length-100mm Diameter-(1x5mm trocar+2x10mm cannulas+1x10mm obturator) + smoke filter
3	IIVT-100-4-B	Optical, Length-100mm Diameter-(2x5mm cannulas+1x5mm obturator+1x10mm trocar+1x12mm trocar) + smoke filter

Note: For more product information, please consult our marketing department.



# STARPORT

## DISPOSABLE SINGLE-INCISION PORT



Website QR Code

**Surgaid** 施爱德®  
Surgaid, an excellent aid for the surgeons

*Surgaid starport percutaneous surgical system, the safety you can feel!*


### Product features



It is widely used in different laparoscopic surgeries such as NOTES, TaTME etc.



**IIIA:Optical**  
easy to observe operation



**IVA:Optical and Rotatable**



**IVC: Rotatable**  
provides flexible operation

### Specifications

No.	Ref No.	Description
1	IIIA-4E-50x150	TPU,4 holes(5/5/12/12mm),diameter-50mm, length-150mm
2	IVA-4E-50x150	TPU,4 holes(5/5/12/12mm),diameter-50mm, length-150mm
3	IVA-4E-70x150	TPU,4 holes(5/5/12/12mm),diameter-70mm, length-150mm
4	IVC-3E-50x70	Silicon, 3 holes(5/12/12mm),diameter-50mm,length-70mm
5	IVC-4D-70x70	Silicon, 4 holes(5/5/10/12mm),diameter-70mm, length-70mm

Note: For more product information, please consult our marketing department.



# STARPORT

## DISPOSABLE SINGLE-INCISION PORT (FLAT TYPE)




Website QR Code

**Surgaid** 施爱德®  
Surgaid, an excellent aid for the surgeons

*Avoiding the spread of tumor cells around the surgical incision, the safety you can feel!*

### Product features



#### Flexible adaption with good sealability

Collocating with Surgaid common trocar or pediatric trocar to supply an operative platform. Meanwhile, the number of holes is selected according to surgical needs to ensure good sealability.

#### Two kinds of models for different choices

TPU material gives the benefits of adjustable height based on different surgery.

Silicone material with fixed height provide effective support.

#### Safe and minimally invasive to protect the incision

Single-incision surgery reduces the number of surgical puncture.

### Specifications

No.	Ref No.	Description
1	IA-3A-50x150	TPU,3 holes(5/5/10mm),diameter-50mm, length-150mm
2	IC-3A-25x30	Silicon,3 holes(5/5/10mm),diameter-25mm, length-30mm
3	IC-3A-50x70	Silicon,3 holes(5/5/10mm),diameter-50mm, length-70mm

Note: For more product information, please consult our marketing department.