

SURGAID

Intelligent Insufflator system

TORNADO™ INSUFFLATOR



【微信公众号】



【外科视界小程序】

Surgaid施爱德®
外科医生好帮手

Edition 202305

Intelligent constant pressure Safe and valve-free

■ 产品特点 Product Features

- Constant pneumoperitoneum pressure, intelligent regulation, meeting the needs of endoscopic and robotic surgery
- Automatic smoke exhaust, mist removal and filtration, maintain a clear surgical field, and create a smoke-free operating room.
- Valveless design, barrier-free entry and exit of instruments and specimens, improving surgical efficiency



TORNADO Mode (Valveless airtight Mode)




- Patented product, realizing valveless pneumoperitoneum and intelligent maintenance of pneumoperitoneum.
- Used together with the visual valveless trocar, the abdominal cavity access without mechanical valve is realized, the instrument is unobstructed, the frequency of cleaning the lens is reduced, and the specimen can be easily taken out.
- Visual operation, expand the surgical field of view, improve surgical safety.

Smoke Extraction Mode



- Circulation gas injection mode, automatic continuous smoke removal.
- The micron filter screen filters all kinds of ultrafine particles, improving the safety of surgery.



Standard Mode

- Used with conventional trocars, high-flow ventilation, stable pneumoperitoneum, to meet the needs of clinical
- Real-time monitoring of abdominal pressure, 1-2 seconds/time, less fluctuation of abdominal pressure, more stable pneumoperitoneum

■ Specifications

NO.	REF NO.	Description
1	QFG-L2000	Insufflation tube, Length-2000mm
2	QFG-L2500	Insufflation tube, Length-2500mm
3	QFG-L3000	Insufflation tube, Length-3000mm
4	QF-A	Insufflation device

Note: For more product information, please consult our marketing department.