


SURGAID

REUSABLE LAPAROSCOPIC INSTRUMENTS (CLIP APPLIER AND REMOVER)



Website QR Code

www.surgaid-medical.com

Surgaid[®]

Surgaid, an excellent aid for the surgeons

Edition 202512

Applier more powerful , remover less effort !

■ Product Introduce

Appliers and clip removers are laparoscopic surgical instruments. They are composed of components such as the jaw, shaft, rotating wheel, and handle.











The product is made from stainless steel and is reusable.

As a delivery device, it is used to apply and remove ligation clips through a trocar of a specific size.

■ Product Features

- ✓ Special clamp head and clamp handle design, more safe and easy to clamp and remove ligature clips;
- ✓ 360°rotating wheel design, no dead angle for instrument use;
- ✓ 10mm working diameter, instruments in and out of the cannula without obstacles;
- ✓ 330, 400mm working length, suitable for most luminal surgery patients;
- ✓ Re-use of instruments, more economical cost, more resource-saving.

■ Product Model

Name	Model	Description	Adaptive Ligature Clip		PCS/BOX
CLIP APPLIERS	F1A01-5x330 F1A01-5x400		S		1 pcs/box
	F1A02-10x330 F1A02-10x400		M		1 pcs/box
	F1A03-10x330 F1A03-10x400		L		1 pcs/box
CLIP REMOVERS	F3A01-5x330 F3A01-5x400		S		1 pcs/box
	F3A01-10x330 F3A01-10x400		M L		1 pcs/box

Note: 5, 10 for the working diameter of the instrument 5, 10mm;330, 400 are the working lengths of the instruments 330, 400 mm.

SURGAID MEDICAL(XIAMEN)CO.,LTD.

No. 18, Shanbianhong East Road, Haicang District, Xiamen 361026, P.R.China
E-mail: info@surgaid.com

Website: www.surgaid-medical.com

Tel: 0592-6885029

Fax: 0592-6587078