




## DISPOSABLE SMOKE FILTER




Website QR Code

**Surgaid** 施爱德®  
Surgaid, an excellent aid for the surgeons

*Breathe healthy and enjoy safety all the time*

### Effective Filter

Contains activated carbon to absorb harmful gases.




### Universal Connector

High compatibility to connect to the trocar.




### Moisture Management

Built-in fluid trap to manage moisture throughout entire procedure.



### Adjustable flow rate

0: Closing velocity    1: 50% velocity  
2: 90% velocity    3: 100% velocity



### Specifications

No.	Ref No.	Description
1	T	Smoke Filter
2	IIVT-100-3-A	Optical 1x5mm 2x10mm smoke filter
3	IIVT-100-4-B	Optical 2x5mm 1x10mm 1x12mm smoke filter

Note: For more product information, please consult our marketing department.