



DISPOSABLE LAPAROSCOPIC LINEAR STAPLER I



Website QR Code

Surgaid 施爱德®

Surgaid, an excellent aid for the surgeons

Disposable laparoscopic linear stapler, the safety you can feel!

Product features

The special "droplets" shaping groove design can ensure B-shape stapler staple formation and provide better hemostatic effect and air tightness effect for clinical operation.

Indicate the direction of the blade walking to ensure safe operation without any accident.

Reset the button manually for firing.

Adaptive joint head can swing within 90°. The joint moves flexible and can be operated by single hand conveniently.

3+1 firing model, manual firing back after 3 step of firing, safer and more stable.

Specifications

Product Code	Device Length	Staple Length	Staple Height
WQD	S - Short (270mm) E - Standard (330mm) L - Long (430mm)	45、60、75	A - 2.5 Y - 3.8 C - 3.5 G - 4.1

No.	REF NO.	Device Length	No.	REF NO.	Device Length	No.	REF NO.	Device Length
1	WQDbS-60A	Short	5	WQDbE-60A	Standard	9	WQDbL-60A	Long
2	WQDbS-60C		6	WQDbE-60C		10	WQDbL-60C	
3	WQDbS-60Y		7	WQDbE-60Y		11	WQDbL-60Y	
4	WQDbS-60G		8	WQDbE-60G		12	WQDbL-60G	

Reload I	Staple Length (mm)	Staple Heights
	45、60、75	A - 2.5 (White) C - 3.5 (Blue) Y - 3.8 (Gold) G - 4.1 (Green)

No.	REF NO.	Staple Length	No.	REF NO.	Staple Length	No.	REF NO.	Staple Length
1	ZQ Nb-45A	45	5	ZQ Nb-60A	60	9	ZQ Nb-75A	75
2	ZQ Nb-45C		6	ZQ Nb-60C		10	ZQ Nb-75C	
3	ZQ Nb-45Y		7	ZQ Nb-60Y		11	ZQ Nb-75Y	
4	ZQ Nb-45G		8	ZQ Nb-60G		12	ZQ Nb-75G	

Note: For more product information, please consult our marketing department.