



DISPOSABLE TROCAR KIT



Website QR Code

Surgaid 施爱德®


Surgaid, an excellent aid for the surgeons

Disposable laparoscopic trocar kit, the safety you can feel!

Product features




Appropriate Sealing
offers good air tightness.



Threaded Design control
insertion speed and depth.



Visualization tip allows
excellent identification.



Smooth conical aiming
facilitates easier insertion.



Simple seal detachment
for versatility during specimen
removal.



Fixed button with advanced
fixation keep scope from moving
during procedures.

Specifications

No.	Ref No.	Description
1	NPVS-100-1-5	Optical, 5*100mm
2	NPVS-100-1-10	Optical, 10*100mm
3	NPVS-100-1-12	Optical, 12*100mm
4	NPVS-100-2-A	Optical 2x5mm
5	NPVS-100-2-B	Optical 2x10mm
6	NPVS-100-2-C	Optical 2x12mm
7	NPVS-100-3-A	Optical 1x5mm 2x10mm
8	NPVS-100-4-A	Optical 2x5mm 2x10mm
9	NPVS-100-4-B	Optical 2x5mm 1x10mm 1x12mm
10	NPVS-100-4-C	Optical 2x5mm 2x12mm

Note: For more product information, please consult our marketing department.