


MINILAP-S

DISPOSABLE FASCIA CLOSURE DEVICE



Website QR Code



Surgaid 施爱德®

Surgaid, an excellent aid for the surgeons


施爱德小儿带穿刺微抓钳，您能感觉的安全!

Surgaid minilap-s, the safety you can feel!

专利设计 PATENT DESIGN



1.6mm diameter



Unique Jaw design



Firmly-hold and uninterrupted


Minimally invasive:

1.6mm minimally invasive diameter reduces one surgical incision and one cannula. Prevent additional incision for children with bilateral hernias.

Reduce operation time:

The unique jaw design avoid withdrawing completely and re-puncture that reduces trauma and operation time.

注水分离 SEPARATION OF WATER INJECTION



Safety:

Unique water injection separation avoids the damage of vas deferens.

Note: Specification: minilap-s-100x1.6. For more product information, Please consult our marking department.